

Product datasheet for **HP215774**

USP42 Human qPCR Primer Pair (NM_032172)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84132
Forward Sequence:	ACAGCGACAGTGACCCGAAAGA
Reverse Sequence:	GACAGCAAAGGAGCAGGCATCA
ACCN:	NM_032172 , NM_032172.1 , NM_032172.2 , BC060846 , BC132862 , BC136656 , BF725643 , BM980885 , BM995466 , BU187025 , BU567897 , NM_032172.3
UniProt ID:	Q9H9J4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)