

Product datasheet for **HP215775**

ZNRF3 Human qPCR Primer Pair (NM_032173)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84133
Forward Sequence:	TCCTGGCTTTCTTCGTCGTGGT
Reverse Sequence:	TGCTCTTGAGTTGAACTTCTGG
ACCN:	NM_032173 , NM_032173.2 , NM_032173.3 , NM_032173.1 , BC009252 , BC021570 , BC042614 , BC051344 , BC069019 , BC094857 , BC167805
UniProt ID:	Q9ULT6
Synonyms:	BK747E2.3; RNF203
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)