

## Product datasheet for **HP215797**

### **ANKRD17 Human qPCR Primer Pair (NM\_032217)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	26057
<b>Forward Sequence:</b>	CAAATGGTGGACACCTCGATGTG
<b>Reverse Sequence:</b>	CTAAGTAGCGCACCTTCAC
<b>ACCN:</b>	<u><a href="#">BC004173</a></u> , <u><a href="#">NM_032217</a></u> , <u><a href="#">NM_032217.1</a></u> , <u><a href="#">NM_032217.2</a></u> , <u><a href="#">NM_032217.3</a></u> , <u><a href="#">NM_032217.4</a></u> , <u><a href="#">BC004891</a></u> , <u><a href="#">BC007747</a></u> , <u><a href="#">BC009043</a></u> , <u><a href="#">BC019963</a></u> , <u><a href="#">BC029935</a></u> , <u><a href="#">BC043394</a></u> , <u><a href="#">BC146382</a></u> , <u><a href="#">BM755557</a></u> , <u><a href="#">BM971664</a></u> , <u><a href="#">BU608983</a></u> , <u><a href="#">BU687888</a></u> , <u><a href="#">BX648035</a></u>
<b>UniProt ID:</b>	<u><a href="#">O75179</a></u>
<b>Synonyms:</b>	GTAR; MASK2; NY-BR-16
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)