

Product datasheet for **HP215799**

CHD6 Human qPCR Primer Pair (NM_032221)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84181
Forward Sequence:	GGACCTCATCAACAAGGAAGCC
Reverse Sequence:	TGATGCGGAAGTGTGTCCAGTC
ACCN:	NM_032221 , NM_032221.1 , NM_032221.2 , NM_032221.3 , NM_032221.4 , BC012516 , BC021907 , BC039860 , BC040016 , BC172397 , BF063641 , BM973992 , NM_032221.5
UniProt ID:	Q8TD26
Synonyms:	CHD-6; CHD5; RIGB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)