

Product datasheet for **HP215808**

SGK196 (POMK) Human qPCR Primer Pair (NM_032237)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84197
Forward Sequence:	AACTGCTCACCTTGCTGTCTT
Reverse Sequence:	TGTGCTCCTTCCACTCAGACAG
ACCN:	NM_032237 , NM_032237.1 , NM_032237.2 , NM_032237.3 , NM_032237.4 , BC113703 , BC101548 , BQ365114 , NM_032237.5
UniProt ID:	Q9H5K3
Synonyms:	MDDGA12; MDDGC12; SGK196
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)