

Product datasheet for **HP215839**

RHBDD1 Human qPCR Primer Pair (NM_032276)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84236
Forward Sequence:	GTTGGTTACCCAGGACGGCAAT
Reverse Sequence:	CTTCACTCAGTCCTGCTGTGTAC
ACCN:	<u>NM_032276</u> , <u>NM_032276.2</u> , <u>NM_032276.3</u> , <u>NM_032276.4</u> , <u>NM_032276.1</u> , <u>BC111056</u> , <u>BC111056.1</u> , <u>BC015553</u> , <u>BC027900</u> , <u>BC062636</u> , <u>BC071679</u> , <u>BC089404</u> , <u>BC101262</u> , <u>BC101263</u> , <u>BC101264</u> , <u>BC101265</u> , <u>NM_032276.5</u>
UniProt ID:	<u>Q8TEB9</u>
Synonyms:	RHBDL4; RRP4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)