

Product datasheet for **HP215863**

CHCHD5 Human qPCR Primer Pair (NM_032309)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84269
Forward Sequence:	CAGCGGGACTGTCACTACCTTA
Reverse Sequence:	TGCCTCGTTCTGTGAAGACAC
ACCN:	<u>NM_032309</u> , <u>NM_032309.1</u> , <u>NM_032309.2</u> , <u>NM_032309.3</u> , <u>BC004498</u> , <u>BC004498.2</u> , <u>BM563093</u> , <u>BX115079</u> , <u>NM_032309.4</u>
UniProt ID:	<u>Q9BSY4</u>
Synonyms:	C2orf9; CHTM1; MIC14; MIX14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)