

Product datasheet for **HP215904**

MON1A Human qPCR Primer Pair (NM_032355)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84315
Forward Sequence:	GGAGGCACTTTCCAGCACTATG
Reverse Sequence:	CGCACGAATACTACCTTGTAGCC
ACCN:	NM_032355 , NM_032355.1 , NM_032355.2 , NM_032355.3 , BC006299 , BC006299.2 , BC009459 , BC047022 , BF514313 , BQ946394
UniProt ID:	Q86VX9
Synonyms:	SAND1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)