

## Product datasheet for **HP215910**

### THOC3 Human qPCR Primer Pair (NM\_032361)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84321
Forward Sequence:	TTCGCATCTGGGATGTGAGGAC
Reverse Sequence:	GTTGCCTACAGCAATGGTCTGC
ACCN:	<a href="#">NM_032361</a> , <a href="#">NM_032361.1</a> , <a href="#">NM_032361.2</a> , <a href="#">NM_032361.3</a> , <a href="#">BC006849</a> , <a href="#">BC006849.1</a> , <a href="#">BC066325</a> , <a href="#">BC068499</a> , <a href="#">BG502664</a> , <a href="#">BU634063</a> , <a href="#">BU781626</a>
UniProt ID:	<a href="#">Q96J01</a>
Synonyms:	hTRES45; THO3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)