

## Product datasheet for **HP215912**

### DNAJC14 Human qPCR Primer Pair (NM\_032364)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	85406
Forward Sequence:	GCCTTCAAGGTTTTGCGAGCAG
Reverse Sequence:	TGCCTCCTTGAGGTCATCTTGC
ACCN:	<a href="#">BC002534</a> , <a href="#">BC016941</a> , <a href="#">BD025734</a> , <a href="#">NM_032364</a> , <a href="#">NM_032364.3</a> , <a href="#">NM_032364.4</a> , <a href="#">NM_032364.5</a> , <a href="#">BC117146</a> , <a href="#">BC080655</a> , <a href="#">BC117148</a>
UniProt ID:	<a href="#">Q6Y2X3</a>
Synonyms:	DNAJ; DRIP78; HDJ3; LIP6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)