

Product datasheet for **HP215925**

HPS3 Human qPCR Primer Pair (NM_032383)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84343
Forward Sequence:	TGTAGCTGAGGCTGGGTGGAAT
Reverse Sequence:	CTGCACTAACAGATGAGCCTC
ACCN:	NM_032383 , NM_032383.1 , NM_032383.2 , NM_032383.3 , NM_032383.4 , BC040359 , BC040359.1 , BC016901 , BC022062 , BG618578 , NM_032383.5
UniProt ID:	Q969F9
Synonyms:	BLOC2S1; SUTAL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)