

## Product datasheet for **HP215953**

### SLX4 Human qPCR Primer Pair (NM\_032444)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84464
Forward Sequence:	AACCCACACCTGAGTGATGTCC
Reverse Sequence:	TCCTCTACAGCGGAGAAGCCTT
ACCN:	<a href="#">NM_032444</a> , <a href="#">NM_032444.1</a> , <a href="#">NM_032444.2</a> , <a href="#">NM_032444.3</a> , <a href="#">BC029952</a> , <a href="#">BC036335</a> , <a href="#">BC104222</a> , <a href="#">BC104223</a> , <a href="#">BC112989</a> , <a href="#">BC112990</a> , <a href="#">BE887716</a> , <a href="#">BM541728</a> , <a href="#">NM_032444.4</a>
UniProt ID:	<a href="#">Q8IY92</a>
Synonyms:	BTBD12; FANCP; MUS312
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)