

Product datasheet for **HP215957**

FAM120B Human qPCR Primer Pair (NM_032448)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84498
Forward Sequence:	TGCCTAGCTGTCAAGGAGTGGT
Reverse Sequence:	GCCGAACCTGTATTTCTGGCTC
ACCN:	<u>NM_032448</u> , <u>NM_032448.1</u> , <u>NM_032448.2</u> , <u>BC012177</u> , <u>BC012177.1</u>
UniProt ID:	<u>Q96EK7</u>
Synonyms:	CCPG; dj894D12.1; KIAA1838; PGCC1; SAN1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)