

Product datasheet for **HP215971**

MCM8 Human qPCR Primer Pair (NM_032485)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84515
Forward Sequence:	CTGTGTGTCGAGGCAGGTCATT
Reverse Sequence:	TCGTGGAATCCGACCTGCTTCT
ACCN:	<u>NM_032485</u> , <u>NM_032485.1</u> , <u>NM_032485.2</u> , <u>NM_032485.3</u> , <u>NM_032485.4</u> , <u>NM_032485.5</u> , <u>BC080656</u> , <u>BC080656.1</u> , <u>BC005170</u> , <u>BC008830</u> , <u>BC101054</u> , <u>BC101055</u> , <u>BC101056</u> , <u>BC101057</u> , <u>BM664368</u> , <u>NM_032485.6</u>
UniProt ID:	<u>Q9UJA3</u>
Synonyms:	C20orf154; dJ967N21.5; POF10
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)