

Product datasheet for **HP215978**

AP1M1 Human qPCR Primer Pair (NM_032493)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8907
Forward Sequence:	GAATGACCGCACCATCTCCTTC
Reverse Sequence:	GAGTGCTTCTCGATCACCGACT
ACCN:	<u>NM_032493</u> , <u>NM_032493.2</u> , <u>NM_032493.3</u> , <u>BC017469</u> , <u>BC017469.2</u> , <u>BQ014582</u> , <u>BU729646</u> , <u>NM_032493.4</u>
UniProt ID:	<u>Q9BXS5</u>
Synonyms:	AP47; CLAPM2; CLTNM; MU-1A; mu1A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)