

Product datasheet for **HP216007**

KIRREL 3 (KIRREL3) Human qPCR Primer Pair (NM_032531)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84623
Forward Sequence:	CCTCGTCCTTATGGCAACCATC
Reverse Sequence:	AGAGGCTGGTTCCTGTGGACA
ACCN:	<u>NM_032531</u> , <u>NM_032531.1</u> , <u>NM_032531.2</u> , <u>NM_032531.3</u> , <u>BC101801</u> , <u>BC101801.1</u> , <u>BC042432</u> , <u>BC065255</u> , <u>BC071838</u> , <u>BC101775</u> , <u>BM671831</u> , <u>NM_032531.4</u>
UniProt ID:	<u>Q8IZU9</u>
Synonyms:	KIRRE; MRD4; NEPH2; PRO4502
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

