

Product datasheet for **HP216015**

SCOC Human qPCR Primer Pair (NM_032547)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	60592
Forward Sequence:	GCTGAAAATCAAGTGGAACTGGAG
Reverse Sequence:	CTTCCTTAAGTGCATCTACTCTTG
ACCN:	NM_032547 , NM_032547.1 , NM_032547.2 , BC062684 , BC016511 , BG151522 , BU070009 , NM_032547.3
UniProt ID:	Q9UIL1
Synonyms:	HRIHFB2072; SCOCO; UNC-69
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)