

Product datasheet for **HP216019**

DAB2IP Human qPCR Primer Pair (NM_032552)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	153090
Forward Sequence:	TCATCGCCAAGGTCACCCAGAA
Reverse Sequence:	CGCTGCATGTTGGTCCACTCAT
ACCN:	NM_032552 , NM_032552.1 , NM_032552.2 , NM_032552.3 , BC011496 , BC037314 , BC040911 , BC146762 , BM545676 , BM546615
UniProt ID:	Q5VWQ8
Synonyms:	AF9Q34; AIP-1; AIP1; DIP1/2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)