

Product datasheet for **HP216065**

CAPS2 Human qPCR Primer Pair (NM_032606)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84698
Forward Sequence:	GGAAGCCTTGTTGACTCTGATGA
Reverse Sequence:	ATGGTGTGACCCAGCCAAGAGT
ACCN:	NM_032606 , NM_032606.1 , NM_032606.2 , NM_032606.3 , BC033860 , BC056672 , BC131724 , BG186487 , NM_032606.5
UniProt ID:	Q9BXY5
Synonyms:	UG0636c06
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)