

Product datasheet for **HP216070**

BEX2 Human qPCR Primer Pair (NM_032621)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84707
Forward Sequence:	TTGGAGAGCCACAGGCAAGGAT
Reverse Sequence:	AGTGCTGACTGCCCGCAAATA
ACCN:	NM_032621 , NM_032621.1 , NM_032621.2 , NM_032621.3 , BC015522 , BC015522.1 , BC050651 , BI670202 , BU569325 , BU596902 , NM_032621.4
UniProt ID:	Q9BXY8
Synonyms:	BEX1; DJ79P11.1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)