

Product datasheet for **HP216074**

PAPOLA Human qPCR Primer Pair (NM_032632)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10914
Forward Sequence:	GTGCTACGAAGACCAGTCCATTG
Reverse Sequence:	GTGGCACAGAACTTCAGTAGCC
ACCN:	<u>BC015150</u> , <u>NM_032632</u> , <u>NM_032632.1</u> , <u>NM_032632.2</u> , <u>NM_032632.3</u> , <u>NM_032632.4</u> , <u>BC036014</u> , <u>BC036014.1</u> , <u>BC000927</u> , <u>BX089817</u> , <u>BX161482</u> , <u>BX248301</u> , <u>BX248753</u>
UniProt ID:	<u>P51003</u>
Synonyms:	PAP; PAP-alpha
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)