

Product datasheet for **HP216083**

CNDP1 Human qPCR Primer Pair (NM_032649)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84735
Forward Sequence:	CCTGTGGATCAGCCAAAGGAAG
Reverse Sequence:	GATGCCACCAAAGTTCCTGAG
ACCN:	NM_032649 , NM_032649.1 , NM_032649.2 , NM_032649.3 , NM_032649.4 , NM_032649.5 , BC110295 , BC110295.1 , BC004271 , BC113512 , BC117122 , BG202268 , BI597577 , BM662490 , NM_032649.6
UniProt ID:	Q96KN2
Synonyms:	CN1; CPGL2; HsT2308
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)