

## Product datasheet for **HP216109**

### KMT5C Human qPCR Primer Pair (NM\_032701)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84787
Forward Sequence:	ACGAGGTGACATGCTTCTACGG
Reverse Sequence:	CTGGTTCGAAAGCTCCTTCAC
ACCN:	<a href="#">BC005842</a> , <a href="#">NM_032701</a> , <a href="#">NM_032701.1</a> , <a href="#">NM_032701.3</a> , <a href="#">BC005842.2</a> , <a href="#">BC044889</a> , <a href="#">BC044889.1</a> , <a href="#">BC019313</a> , <a href="#">NM_032701.4</a>
UniProt ID:	<a href="#">Q86Y97</a>
Synonyms:	Suv4-20h2; SUV420H2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)