

## Product datasheet for **HP216117**

### MFSD9 Human qPCR Primer Pair (NM\_032718)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 84804  |
| Forward Sequence: | TTCTCGGAGCAGCCACCAATGT   |
| Reverse Sequence: | GAAGTGTCGGATTACAAGCGGC   |
| ACCN:             | <a href="#">NM_032718</a> , <a href="#">NM_032718.1</a> , <a href="#">NM_032718.2</a> , <a href="#">NM_032718.3</a> , <a href="#">NM_032718.4</a> , <a href="#">BC006242</a> ,<br><a href="#">BC006242.2</a> , <a href="#">BC009310</a>  |
| UniProt ID:       | <a href="#">Q8NBP5</a>   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)