

## Product datasheet for **HP216156**

### TIM 3 (HAVCR2) Human qPCR Primer Pair (NM\_032782)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84868
Forward Sequence:	GACTCTAGCAGACAGTGGGATC
Reverse Sequence:	GGTGGTAAGCATCCTTGGAAAGG
ACCN:	<a href="#">BC020843</a> , <a href="#">NM_032782</a> , <a href="#">NM_032782.1</a> , <a href="#">NM_032782.2</a> , <a href="#">NM_032782.3</a> , <a href="#">NM_032782.4</a> , <a href="#">BC020843.1</a> , <a href="#">BC063431</a> , <a href="#">BC063431.1</a> , <a href="#">BP291914</a>
UniProt ID:	<a href="#">Q8TDQ0</a>
Synonyms:	CD366; HAVcr-2; KIM-3; SPTCL; Tim-3; TIM3; TIMD-3; TIMD3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)