

Product datasheet for **HP216170**

ZDH12 (ZDHHC12) Human qPCR Primer Pair (NM_032799)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84885
Forward Sequence:	GAATCACGCTGGTGCTCTTCCT
Reverse Sequence:	AGTGACACAGCGAGGTAGAGCA
ACCN:	BC048251 , NM_032799 , NM_032799.1 , NM_032799.2 , NM_032799.3 , NM_032799.4 , BC048251.1 , BC009280 , BC009280.2 , BG767826 , BM043458 , BM745076 , BU506687
UniProt ID:	Q96GR4
Synonyms:	DHHC-12; ZNF400
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)