

Product datasheet for **HP216207**

RITA (RITA1) Human qPCR Primer Pair (NM_032848)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84934
Forward Sequence:	AACCAGAGCGTGGCAAGGA
Reverse Sequence:	GCGACTCATCACAGTAAGACGG
ACCN:	<u>NM_032848</u> , <u>NM_032848.1</u> , <u>NM_032848.2</u> , <u>BC022092</u> , <u>BC022092.2</u> , <u>BP212183</u> , <u>BX366456</u> , <u>NM_032848.3</u>
UniProt ID:	<u>Q96K30</u>
Synonyms:	C12orf52; RITA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)