

Product datasheet for **HP216210**

PWWP3A Human qPCR Primer Pair (NM_032853)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84939
Forward Sequence:	TGTGTCTCCCTCATCACCGACT
Reverse Sequence:	GCTGGATGGATTTTCGGACAGG
ACCN:	<u>BD127600</u> , <u>NM_032853</u> , <u>NM_032853.1</u> , <u>NM_032853.2</u> , <u>NM_032853.3</u> , <u>BC110874</u> , <u>BC110874.1</u> , <u>BC001806</u> , <u>BC008098</u> , <u>BC019585</u> , <u>BC021016</u> , <u>BC082987</u> , <u>BC130443</u> , <u>BC144138</u> , <u>BC144147</u> , <u>BG723598</u> , <u>BI553935</u> , <u>BQ071619</u> , <u>BQ895780</u> , <u>BU541563</u> , <u>BU741593</u> , <u>NM_032853.5</u>
UniProt ID:	<u>Q2TAK8</u>
Synonyms:	EXPAND1; HSPC211; MUM-1; MUM1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)