

Product datasheet for **HP216214**

LACTB Human qPCR Primer Pair (NM_032857)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	114294
Forward Sequence:	GTGCAGGAAGAAAACGAGCCAG
Reverse Sequence:	CAGAAGGTCACCCACTGTAGAC
ACCN:	<u>NM_032857</u> , <u>NM_032857.1</u> , <u>NM_032857.2</u> , <u>NM_032857.3</u> , <u>NM_032857.4</u> , <u>BC067288</u> , <u>BC037198</u> , <u>BG250900</u> , <u>BU624474</u> , <u>NM_032857.5</u>
UniProt ID:	<u>P83111</u>
Synonyms:	G24; MRPL56
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)