

Product datasheet for **HP216216**

ABHD13 Human qPCR Primer Pair (NM_032859)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84945
Forward Sequence:	GTCTTGGGCTCTCTGCCGTATT
Reverse Sequence:	CGTGAAGAGGATGGCTGTTCTG
ACCN:	NM_032859 , NM_032859.1 , NM_032859.2 , BC022566 , BC022566.1 , BC018608 , BC070226 , BQ424339 , BX647428 , NM_032859.3
UniProt ID:	Q7L211
Synonyms:	bA153I24.2; BEM46L1; C13orf6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)