

Product datasheet for **HP216259**

RAB11FIP4 Human qPCR Primer Pair (NM_032932)

Product data:

| | |
|--------------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 84440 |
| Forward Sequence: | CCTGCTCAATGACTTGGAAGCC |
| Reverse Sequence: | ACAGGTCTTCGGTGTGCCATT |
| ACCN: | BC006240 , BC006563 , NM_032932 , NM_032932.2 , NM_032932.3 , NM_032932.4 , NM_032932.5 , BC093914 , BC093914.1 , BC004515 , BC023618 , BC047674 , BC094803 , BC101517 , BC131530 , BI461640 |
| UniProt ID: | Q86YS3 |
| Synonyms: | FIP4-Rab11; RAB11-FIP4 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)