

Product datasheet for **HP216262**

TMEM60 Human qPCR Primer Pair (NM_032936)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	85025
Forward Sequence:	AGCACCTTGGAAGTGGTTCCTC
Reverse Sequence:	GGTCAAAGCCAGACTTACACCG
ACCN:	<u>NM_032936</u> , <u>NM_032936.1</u> , <u>NM_032936.2</u> , <u>NM_032936.3</u> , <u>BC065930</u> , <u>BC065930.1</u> , <u>BC012192</u> , <u>BX117806</u> , <u>NM_032936.4</u>
UniProt ID:	<u>Q9H2L4</u>
Synonyms:	C7orf35; DC32
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)