

Product datasheet for **HP216290**

OR2B2 Human qPCR Primer Pair (NM_033057)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	81697
Forward Sequence:	AGCTGCATCCTGGATTAGTGCC
Reverse Sequence:	CACAGGACAACTTGAGCAGAGC
ACCN:	<u>NM_033057</u> , <u>NM_033057.1</u> , <u>BC136881</u> , <u>BC136901</u>
UniProt ID:	<u>Q9GZK3</u>
Synonyms:	dj193B12.4; hs6M1-10; OR2B2Q; OR2B9; OR6-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)