

Product datasheet for **HP216306**

FAM40A (STRIP1) Human qPCR Primer Pair (NM_033088)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	85369
Forward Sequence:	GGATGTAACCGCCACAAAGAGG
Reverse Sequence:	TGCAGTTGGCAAACACCAGGTG
ACCN:	NM_033088 , NM_033088.1 , NM_033088.2 , NM_033088.3 , BC094786 , BC094786.1 , BC019919 , BC032644 , BC040142 , BC067239 , BC119814 , BC121793 , NM_033088.4
UniProt ID:	Q5VSL9
Synonyms:	FAM40A; FAR11A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)