

Product datasheet for **HP216312**

STON2 Human qPCR Primer Pair (NM_033104)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	85439
Forward Sequence:	GTGTGGAGGATAAACCGACTGC
Reverse Sequence:	CGTGATTGGCAAATCTGGAAGGC
ACCN:	<u>NM_033104</u> , <u>NM_033104.1</u> , <u>NM_033104.2</u> , <u>NM_033104.3</u> , <u>BC069359</u> , <u>BC069389</u> , <u>BC069502</u> , <u>BC069587</u> , <u>BC117493</u> , <u>BC117495</u> , <u>BX648249</u>
UniProt ID:	<u>Q8WXE9</u>
Synonyms:	STN2; STNB; STNB2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)