

Product datasheet for **HP216321**

NEK9 Human qPCR Primer Pair (NM_033116)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	91754
Forward Sequence:	GGAGACAAAGCCTCCTATCGAC
Reverse Sequence:	ATCCGAAGGCATAGAGCTGACC
ACCN:	<u>NM_033116</u> , <u>NM_033116.1</u> , <u>NM_033116.3</u> , <u>NM_033116.5</u> , <u>BC112101</u> , <u>BC112101.1</u> , <u>BC001714</u> , <u>BC009336</u> , <u>BC093881</u> , <u>BM976692</u> , <u>NM_033116.6</u>
UniProt ID:	<u>Q8TDI9</u>
Synonyms:	APUG; LCCS10; NC; NERCC; NERCC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

