

## Product datasheet for **HP216328**

### PLCZ1 Human qPCR Primer Pair (NM\_033123)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	89869
Forward Sequence:	GCAGGCTACTTTTGGAGAGTCC
Reverse Sequence:	TTCATGGGTTTCCTTAAGGTTCC
ACCN:	<u>NM_033123</u> , <u>NM_033123.1</u> , <u>NM_033123.3</u> , <u>BC125067</u> , <u>NM_033123.4</u>
UniProt ID:	<u>Q86YW0</u>
Synonyms:	NYD-SP27; PLC-zeta-1; PLCzeta; SPGF17
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)