

Product datasheet for **HP216348**

KRTAP2-4 Human qPCR Primer Pair (NM_033184)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	85294
Forward Sequence:	TGCTGCCGCCCATCACCTG
Reverse Sequence:	AGACACAGGCTGGCAGCAGGT
ACCN:	<u>NM_033184</u> , <u>NM_033184.1</u> , <u>NM_033184.2</u> , <u>NM_033184.3</u> , <u>BC063625.1</u> , <u>BC101147</u> , <u>BC101148</u> , <u>BC101149</u>
UniProt ID:	<u>P0C7H8</u>
Synonyms:	KAP2.1B; KAP2.3; KAP2.4; KRTAP2-3; KRTAP2.3; KRTAP2.4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)