

## Product datasheet for **HP216370**

### Pur B (PURB) Human qPCR Primer Pair (NM\_033224)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5814
Forward Sequence:	GATGTGAAGCAGAACGCCAAGG
Reverse Sequence:	GCGTAGTGTCTATGAAGTCGCC
ACCN:	<a href="#">NM_033224</a> , <a href="#">NM_033224.1</a> , <a href="#">NM_033224.2</a> , <a href="#">NM_033224.3</a> , <a href="#">NM_033224.4</a> , <a href="#">BC101735</a> , <a href="#">BC101737</a> , <a href="#">BF964815</a> , <a href="#">BG776957</a> , <a href="#">BM015566</a> , <a href="#">BQ890213</a> , <a href="#">BU154650</a> , <a href="#">BU177276</a> , <a href="#">NM_033224.5</a>
UniProt ID:	<a href="#">Q96QR8</a>
Synonyms:	PURBETA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)