

## Product datasheet for **HP216371**

### CSMD1 Human qPCR Primer Pair (NM\_033225)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64478
Forward Sequence:	CGAAGAGGTACAATCCTGTCCC
Reverse Sequence:	CCAGTTCTGCTCCGTGGCAAAA
ACCN:	<a href="#">NM_033225</a> , <a href="#">NM_033225.1</a> , <a href="#">NM_033225.2</a> , <a href="#">NM_033225.3</a> , <a href="#">NM_033225.4</a> , <a href="#">NM_033225.5</a> , <a href="#">BC156304</a> , <a href="#">BP229721</a> , <a href="#">BX107340</a> , <a href="#">NM_033225.6</a>
UniProt ID:	<a href="#">Q96PZ7</a>
Synonyms:	PPP1R24
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)