

Product datasheet for **HP216438**

IL33 Human qPCR Primer Pair (NM_033439)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	90865
Forward Sequence:	GCCTGTCAACAGCAGTCTACTG
Reverse Sequence:	TGTGCTTAGAGAAGCAAGATACTC
ACCN:	<u>NM_033439</u> , <u>NM_033439.1</u> , <u>NM_033439.2</u> , <u>NM_033439.3</u> , <u>BC047085</u> , <u>BC047085.1</u> , <u>BQ878541</u> , <u>NM_033439.4</u>
UniProt ID:	<u>O95760</u>
Synonyms:	C9orf26; DVS27; IL1F11; NF-HEV; NFEHEV
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)