

Product datasheet for **HP216451**

EPT1 (SELENOI) Human qPCR Primer Pair (NM_033505)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	85465
Forward Sequence:	GTGGATACCAATCCAATTTCTCTG
Reverse Sequence:	GAATACGACCAGCAGAAAGCCAG
ACCN:	NM_033505 , NM_033505.1 , NM_033505.2 , NM_033505.3 , BC015430 , BC021229 , BE963795 , BF669898 , BK001426 , BU608897 , NM_033505.4
UniProt ID:	Q9C0D9
Synonyms:	EPT1; SELI; SEPI; SPG81
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)