

Product datasheet for **HP216467**

HLA DP (HLA-DPA1) Human qPCR Primer Pair (NM_033554)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3113
Forward Sequence:	ATCCAGCGTTCCAACCACACTC
Reverse Sequence:	CGTTGAGCACTGGTGGGAAGAA
ACCN:	<u>BC009956</u> , <u>NM_033554</u> , <u>NM_033554.1</u> , <u>NM_033554.2</u> , <u>NM_033554.3</u> , <u>BC009956.2</u> , <u>BI553587</u> , <u>BP338039</u> , <u>BU687966</u>
UniProt ID:	<u>P20036</u>
Synonyms:	DP(W3); DP(W4); DPA1; HLA-DP1A; HLA-DPA; HLA-DPB1; HLADP; HLASB; PLT1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)