

## Product datasheet for **HP216475**

### Pannexin 2 (PANX2) Human qPCR Primer Pair (NM\_052839)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56666
Forward Sequence:	CCTCACCATCCTGAGCCGAAAC
Reverse Sequence:	GACGGTGCCATCCCTCAAACCT
ACCN:	<a href="#">NM_052839</a> , <a href="#">NM_052839.1</a> , <a href="#">NM_052839.2</a> , <a href="#">NM_052839.3</a> , <a href="#">BC101023</a> , <a href="#">BQ576292</a> , <a href="#">BU732466</a> , <a href="#">BU734128</a> , <a href="#">NM_052839.4</a>
UniProt ID:	<a href="#">Q96RD6</a>
Synonyms:	hPANX2; PX2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)