

Product datasheet for **HP216482**

GADD45GIP1 Human qPCR Primer Pair (NM_052850)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	90480
Forward Sequence:	GCCGAAGAACGCGAATGGTACC
Reverse Sequence:	CACAATCATCTGTGGCATCTTGG
ACCN:	<u>NM_052850</u> , <u>NM_052850.1</u> , <u>NM_052850.2</u> , <u>NM_052850.3</u> , <u>BC069200</u> , <u>BC004944</u> , <u>BC013027</u> , <u>BC013039</u> , <u>BI459923</u> , <u>NM_052850.4</u>
UniProt ID:	<u>Q8TAE8</u>
Synonyms:	CKBBP2; CKbetaBP2; CRIF1; MRP-L59; PLINP; PLINP-1; Plinp1; PRG6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)