

Product datasheet for **HP216510**

PGRPL (PGLYRP2) Human qPCR Primer Pair (NM_052890)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	114770
Forward Sequence:	TGGCTGCCAATGCTACCAAGGA
Reverse Sequence:	CGTAGGTGTGATGCACGTACAAG
ACCN:	<u>NM_052890</u> , <u>NM_052890.1</u> , <u>NM_052890.2</u> , <u>NM_052890.3</u> , <u>BC136551</u> , <u>BC136552</u> , <u>BC144238</u> , <u>NM_052890.4</u>
UniProt ID:	<u>Q96PD5</u>
Synonyms:	HMFT0141; PGLYRPL; PGRP-L; PGRPL; tagL; tagL-alpha; tagL-beta; TAGL-like
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)