

Product datasheet for **HP216523**

Eco1 (ESCO1) Human qPCR Primer Pair (NM_052911)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	114799
Forward Sequence:	GCTGAATACCCTGATGGCAGGA
Reverse Sequence:	GGAATAGCACATTAGTGGAGCCT
ACCN:	<u>NM_052911</u> , <u>NM_052911.1</u> , <u>NM_052911.2</u> , <u>BC036943</u> , <u>BC089426</u> , <u>BG492713</u> , <u>BK001617</u> , <u>NM_052911.3</u>
UniProt ID:	<u>Q5FWF5</u>
Synonyms:	A930014I12Rik; CTF; ECO1; EFO1; ESO1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)