

Product datasheet for **HP216528**

SCAND3 (ZBED9) Human qPCR Primer Pair (NM_052923)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	114821
Forward Sequence:	TGCTAAAGCCATCAGGGACACC
Reverse Sequence:	ACTACTGTGTTTCAGAAAAGGCATG
ACCN:	<u>NM_052923</u> , <u>NM_052923.1</u> , <u>BC110834</u> , <u>BC110834.1</u> , <u>BC111741</u> , <u>BC111742</u> , <u>BX647791</u> , <u>NM_052923.2</u>
UniProt ID:	<u>Q6R2W3</u>
Synonyms:	Buster4; dj1186N24.3; SCAND3; ZFP38-L; ZNF305P2; ZNF452
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)