

Product datasheet for **HP216533**

LY108 (SLAMF6) Human qPCR Primer Pair (NM_052931)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	114836
Forward Sequence:	CTGTGTATGCTTCAGTCACTCATT
Reverse Sequence:	CACGACATTGTCAAGGGCAGTTG
ACCN:	<u>NM_052931</u> , <u>NM_052931.1</u> , <u>NM_052931.2</u> , <u>NM_052931.4</u> , <u>BC113893</u> , <u>BC090928</u> , <u>BC114495</u> , <u>NM_052931.5</u>
UniProt ID:	<u>Q96DU3</u>
Synonyms:	CD352; KALI; KAL1b; Ly108; NTB-A; NTBA; SF2000
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)