

Product datasheet for **HP216538**

ATG4A Human qPCR Primer Pair (NM_052936)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	115201
Forward Sequence:	CCAAGCCAGAAGTGACAACCAC
Reverse Sequence:	GACAGACCTTCAAGTTGAGTTCC
ACCN:	NM_052936 , NM_052936.1 , NM_052936.2 , NM_052936.3 , NM_052936.4 , BC041862 , BC061696 , BM541631 , NM_052936.5
UniProt ID:	Q8WYN0
Synonyms:	APG4A; AUTL2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)